

Validation of a real-time RT-PCR based detection system for Yellow Fever Virus specific RNA

Patricia Sequeira², Hussein El Halas¹, Simone Alves Sampaio², Karin Rottengatter¹, Ana Bispo de Filippis², Leonie-S. Hecht¹

¹altona Diagnostics GmbH, Mörkenstraße 12; Hamburg, Germany, ²Instituto Oswaldo Cruz – Fiocruz, Av. Brasil - 4365 Manguinhos, Rio de Janeiro, Brasil

Introduction

Brazil has recently experienced its largest-recorded yellow fever (YF) outbreak in decades, with more than 2000 confirmed cases and over 600 deaths since December 2016. YF is a severe mosquito-borne viral infection in the tropics, responsible for 29000 to 60000 deaths annually in South America and Africa. Yellow fever virus (YFV) is the prototype of the genus *Flavivirus*, which comprises around 70 different arthropod-borne viruses.

No specific treatments have been found to benefit patients with yellow fever, only supportive care to treat dehydration, respiratory failure, and fever. But a YF vaccine is available, which is a live attenuated viral vaccine from the 17D lineage, which elicits a rapid, exceptionally strong, and markedly durable adaptive immune response.

The clinical diagnosis of YF is difficult because of the similarities of the symptoms to those of a wide range of diseases, including dengue fever, other hemorrhagic viral diseases, leptospirosis, viral hepatitis, and malaria; hence laboratory confirmation is essential. Serologic diagnosis is best accomplished using an enzyme-linked immunosorbent assay (ELISA) for IgM. Serological techniques often cross-react among flavivirus infections, thus, the use of real-time reverse transcription polymerase chain reaction (RT-PCR) should be prioritized. Viral RNA can be detected in serum during the first 10 days from symptom onset (viremic phase) using real-time RT-PCR. Samples that test negative by RT-PCR should subsequently be tested by serology.

Here we describe the verification and validation of a new IVD-CE marked real-time RT-PCR assay (RealStar[®] Yellow Fever Virus RT-PCR Kit; altona Diagnostics) for the detection of YFV specific RNA.

Material and Methods

Here we describe the verification and validation of a newly developed yellow fever virus RT-PCR, which comprised of a specific system for the detection of YFV specific RNA (Figure 1 A) and an Internal Control (IC; Figure 1 B) system that monitors the efficiency of nucleic acid extraction process and possible inhibitory effects during PCR.

The RealStar[®] Yellow Fever Virus RT-PCR Kit limit of detection was determined using probit analysis after testing replicates of limited dilutions of quantified *in vitro* transcribed RNA containing the RT-PCR target sequence (Table 1). The analytical sensitivity is 0.69 copies/ μ l [95% confidence interval (CI): 0.41 - 1.56 copies/ μ l], Figure 2.

The analytical specificity of the YFV RT-PCR was evaluated by testing a panel of genomic RNA/DNA extracted from different pathogens that are related to YFV and/or can cause symptoms similar to YF virus (Table 2).

For the diagnostic validation 45 serum samples, pre-characterized with a reference real-time PCR assay, were used at the Flavivirus Laboratory at Fiocruz which is a Regional Reference for the Brazilian MoH.

Aim of the study

The aim of this study was to compare the newly developed yellow fever virus RT-PCR assay with the yellow fever virus RT-PCR assay established at the Flavivirus Laboratory FIOCRUZ/ Brazil (based on Domingo *et al.*, 2012) on a total of 30 serum samples from patients with yellow fever (YF) virus infection, retrospectively. Nucleic acid was extracted and the eluates were tested in parallel using the two different RT-PCRs. In addition, 15 individual serum samples from non-infected persons were tested.

Results and Conclusion

All 30 positive samples for YF virus RNA with the reference were also tested positive using the newly developed YFV RT-PCR assay. Out of the 15 negative samples for YF virus RNA using the reference assay, 15 samples were tested negative with the altona Diagnostics yellow fever virus RT-PCR (Table 3).

Our verification and validation data of the RealStar[®] Yellow Fever Virus RT-PCR Kit indicate that the assay is appropriate for the sensitive and specific detection of YFV specific RNA. It can be a useful tool in patient and epidemiological management in YF endemic regions.

Contact

leonie.hecht@altona-diagnostics.com

altona Diagnostics GmbH, Mörkenstr. 12, 22767 Hamburg, Germany




Figure 1. YFV negative samples (no amplification) and different concentrations of YFV specific RNA, detected in the FAM-Channel (A) and the regarding Internal Control signals (IC) detected in the JOE/VIC-Channel (B).

Table 1: RT-PCR results used for the calculation of the analytical sensitivity with respect to the detection of YFV specific RNA

Input Conc. [copies/ μ l]	Number of Replicates	Number of Positives	Hit Rate [%]
31.600	24	24	100
10.000	24	24	100
3.160	24	24	100
1.000	24	24	100
0.316	24	21	87.5
0.100	24	9	37.5
0.032	24	4	16.7
0.010	24	2	8.3
0.003	24	0	0




Figure 2. The analytic sensitivity is 0.69 copies/ μ l [95% confidence interval (CI): 0.41 - 1.56 copies/ μ l]

Table 2: Analytic specificity results

Pathogens	FAM (Yellow Fever Virus)	JOE/VIC (Internal Control)
Chikungunya Virus	No Ct	Valid
Crimean-Congo Hemorrhagic Fever virus	No Ct	Valid
Dengue Virus Serotype 1	No Ct	Valid
Dengue Virus Serotype 2	No Ct	Valid
Dengue Virus Serotype 3	No Ct	Valid
Dengue Virus Serotype 4	No Ct	Valid
Ebola Virus	No Ct	Valid
Hepatitis C Virus	No Ct	Valid
Japanese encephalitis virus	No Ct	Valid
Lassa Virus	No Ct	Valid
Marburg Virus	No Ct	Valid
Murray Valley encephalitis virus	No Ct	Valid
<i>Plasmodium falciparum</i>	No Ct	Valid
West Nile Virus	No Ct	Valid
Zika Virus	No Ct	Valid
Positive Control	32.25	Valid
NTC	No Ct	Valid

Table 3: Diagnostic validation

Total number of samples: 45	YF Virus real-time RT-PCR (based on Domingo <i>et al.</i> , 2012)	
	POSITIVE	NEGATIVE
RealStar [®] YFV RT-PCR Kit 1.0	POSITIVE	30
	NEGATIVE	15