

Adeno

- Detection and quantification of human Adenovirus (HAdV) specific DNA
- Detection of all known human pathogen Adenovirus species
- Ready-to-use kit including Internal Control and Quantification Standards
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test

RealStar® Adeno PCR Kit 1.0 Prod. No 301013



For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform, CFX96™/Dx Real-Time System

always a drop ahead.



RealStar® Adenovirus PCR Kit



HAdV Quantification Standards QS1 (10^4 copies/ μ l) to QS4 (10^1 copies/ μ l) detected in the FAM detection channel [A]. The Internal Control (IC) is detected in the JOE detection channel [B].

