

Dengue

- Detection of Dengue Virus specific RNA
- Detection of all four Dengue Virus Types 1–4
- Ready-to-use kit including Internal and Positive Control
- Compatible with various Real-Time PCR platforms



RealStar® Dengue RT-PCR Kit 2.0 Prod. No 282003

For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS, Mx 3005P™ QPCR System, VERSANT® kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene® 6000, Rotor-Gene® Q 5/6 plex Platform, CFX96™ Real-Time System

always a drop ahead.

RealStar® Dengue RT-PCR Kit



Dengue Virus specific RNA detected in the FAM detection channel [A]. Red: Dengue Type 1; Blue: Dengue Type 4; Purple: Dengue Type 3; Dark Red: Dengue Type 2; Green: Positive Control; Blue: Negative Control. The Internal Control (IC) is detected in the JOE detection channel [B].

