

## EBV

- Detection and quantification of Epstein-Barr Virus (EBV) specific DNA
- Parallel testing within Herpes- and Polyomavirus test kit panel
- Ready-to-use kit including Internal Control and Quantification Standards
- Quantification in International Units (IU)
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test

### RealStar® EBV PCR Kit 1.0

Prod. No 131013

**For use with:**

ABI Prism® 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform

### RealStar® EBV PCR Kit 1.2



Prod. No 131212

**For use with:**

SmartCycler® II, LightCycler® 1.2 / 1.5 / 2.0



## RealStar® EBV PCR Kit



Dilution series of EBV specific DNA, ranging from  $10^0$  to  $10^7$  copies/ $\mu$ l, detected in the FAM detection channel [A]. The Internal Control is detected in the JOE detection channel [B].

