

HEV



- Detection of Hepatitis E Virus (HEV) specific RNA
- Parallel testing within Hepatitis RT-PCR test kit panel
- Ready-to-use kit including Internal and Postive Control
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test

RealStar® HEV RT-PCR Kit 1.0
Prod. No 271013

For use with:

ABI Prism® 7500 SDS, m2000rt,
Mx 3005P™ QPCR System,
VERSANT™ kPCR Molecular System AD,
LightCycler® 480 Instrument II,
Rotor-Gene™ 3000/6000,
Rotor-Gene™ Q 5/6 plex Platform
CFX96™/Dx Real-Time System

RealStar® HEV RT-PCR Kit



Dilution series of HEV specific RNA, ranging from 10^1 to 10^4 copies/ μ l, detected in the FAM detection channel [A]. The Internal Control (IC) is detected in the JOE detection channel [B].

