


Causative agent of COVID-19



RealStar® SARS-CoV-2 RT-PCR Kit 1.0*

A screen and confirm test, based on real-time PCR technology. Ready-to-use kit including Mastermix, Positive Control and Internal Control.

Screen and confirm in one tube!

Dual target assay to rapidly screen lineage B-betacoronaviruses and confirm the SARS-CoV-2 specific RNA in **only one reaction**.

Target	Signal	Detection
lineage B-betacoronaviruses (E gene)	FAM™	✓
SARS-CoV-2 specific (S gene)	Cy5	✓
Internal Control	JOE™	✓

RealStar® SARS-CoV-2 RT-PCR Kit 1.0*

Prod. No 821003

96 reactions

*For Research Use Only!



Save time and test in parallel!

RealStar® RSV RT-PCR Kit 3.0

Prod. No 193013

96 reactions

Respiratory Syncytial Virus A and B

RealStar® Influenza Screen & Type RT-PCR Kit 4.0

Prod. No 164013

96 reactions

Human Influenza A, B and SwineFlu (H1N1)pdm09